

SPESolid Phase
Extraction

Application Note

World Leaders in Sample Preparation

SPE181**ISOLATION OF PHENOLIC COMPOUNDS
FROM WATER**

Spe-ed™ Cartridge	Cat. No. 4264- <i>Spe-ed Advanta</i> , 100mg/1mL
Cartridge Conditioning	2mL of methanol. 2mLs water.
Sample Addition	Aspirate 20mLs of sample through cartridge.
Cartridge Wash	2mL of water. Air dry cartridge for 1 minute.
Analyte Elution	1mL of ethyl acetate.

Note: *Since sample matrix interferences and concentrations may vary from sample to sample, it may be necessary to adjust the wash and elution solvent/solution strength and/or volume to optimize isolation.*

The logo for Applied Separations, featuring the company name in a bold, sans-serif font with horizontal lines above the word "Applied".

930 Hamilton Street · Allentown, PA 18101
610-770-0900 · 610-740-5520 (fax)
www.appliedseparations.com